



Call to receive our **FREE** Buyer's Guide Catalog

M2250 SERIES

- Monocular head rotates 360°
- Locked-on 10x eyepiece & pointer
- 4x, 10x, & 40x DIN objectives
- Separate coarse & fine focusing
- Stage 4 3/4" x 4 1/2" with stage clips
- Tungsten or LED illumination
- Limited lifetime warranty
- Microscopes by Swift Optical



Catalog #	Description	List Price	Sale
M2251B	Corded Tungsten	\$259 ⁰⁰	\$239⁰⁰
M2251CL	Corded LED	\$275 ⁰⁰	\$249⁰⁰
M2251C	Cordless LED	\$299 ⁰⁰	\$269⁰⁰

ST SERIES

- Paired 10x widefield eyepieces
- Inclined 45° binocular with diopter adjustment
- Parfocalled & parcentered achromatic objectives
- Adjustable post mounted head rotates 180°
- Both top & bottom lighting, rack & pinion focusing
- Interpupillary distance adjustable 51-75 mm
- Microscopes by Walter Products




Catalog #	Description	List Price	Sale
WA-ST1310	10x & 30x / Corded LED	\$259 ⁰⁰	\$209⁰⁰
WA-ST1310RC	10x & 30x / Cordless LED	\$279 ⁰⁰	\$229⁰⁰
WA-ST2410	20x & 40x / Corded LED	\$259 ⁰⁰	\$209⁰⁰
WA-ST2410RC	20x & 40x / Cordless LED	\$279 ⁰⁰	\$229⁰⁰

450 SERIES

NATIONAL

- Fixed 45° inclined head, 12" tall
- Widefield 10x eyepieces, diopter adjustment
- Achromatic objectives parfocalled/parcentered
- LED illumination, 4-position light control: top only, bottom only, top & bottom, or off
- Rack & pinion focusing
- Includes two 80 mm stage plates: frosted & reversible black/white plastic contrast plates
- Microscopes by National Optical



Catalog #	Description	List Price	Sale
NT452TBL-10LED	1x & 3x	\$285 ⁰⁰	\$259⁰⁰
NT453TBL-10LED	2x & 4x	\$285 ⁰⁰	\$259⁰⁰

SM100 SERIES

- W10x eyepiece with diopter
- Top & bottom illumination
- Slip-clutch system
- Stage plates included
- Microscopes by Swift Optical



Catalog #	Description	List Price	Sale
SM101	10x & 30x Corded LED	\$389 ⁰⁰	\$359⁰⁰
SM101C	10x & 30x Cordless LED	\$399 ⁰⁰	\$369⁰⁰
SM105	Zoom 10x-30x Corded LED	\$559 ⁰⁰	\$519⁰⁰
SM105C	Zoom 10x-30x Cordless LED	\$575 ⁰⁰	\$525⁰⁰

M3600 SERIES

- Monocular head
- Pointmaster® eyepiece (W10x)
- Separate coarse & fine focusing
- Locked-on eye tube & stage clips
- Corded or cordless LED lighting
- Reversed nosepiece with DIN objectives
- 100xR objective, mechanical stage & Abbe condenser on select models
- Limited lifetime warranty
- Microscopes by Swift Optical



Catalog #	Description	List Price	Sale
M3601	4x, 10x, 40xR Corded LED	\$379 ⁰⁰	\$349⁰⁰
M3601C	4x, 10x, 40xR Cordless LED	\$389 ⁰⁰	\$359⁰⁰
M3602-4	4x, 10x, 40xR, 100xR Mechanical Stage Abbe, Corded LED	\$495 ⁰⁰	\$445⁰⁰
M3602C-4	4x, 10x, 40xR, 100xR Mechanical Stage Abbe, Cordless LED	\$499 ⁰⁰	\$449⁰⁰

More models available - call us or visit our website for details!

Need help choosing a microscope? Call us for recommendations.

M3700 SERIES

- Advanced binocular with 10x widefield eyepieces, 45° inclined, 360° rotatable head
- Reversed nosepiece with semi-plan parfocalled & parcentered objectives: 4xD, 10xD, 40xRD & 100xRD
- Built-in low-drive mechanical stage
- Coaxial focusing with Abbe condenser
- Cordless LED with dimmer switch
- Microscope by Swift Optical



Catalog #	Description	List Price	Sale
M3702CB-4	4x, 10x, 40xR, 100xR Mechanical Stage Abbe, Cordless LED	\$799 ⁰⁰	\$729⁰⁰

167 SERIES

- Monocular or binocular inclined head rotates 360°
- Widefield 10x eyepiece with pointer
- LED variable illumination
- Coaxial coarse & fine focusing
- Large 140 mm x 135 mm rackless mechanical stage
- DIN achromatic, parfocalled, parcentered, color-coded objectives: 4x, 10x, 40xR & 100xR, reversed nosepiece
- 1.25 Abbe condenser with rack & pinion focusing
- Microscopes by National Optical, limited lifetime warranty

NATIONAL



Catalog #	Description	List Price	Sale
NT167-SP	Monocular	\$626 ⁰⁰	\$559⁰⁰
NT168-SP	Binocular	\$815 ⁰⁰	\$739⁰⁰





2053 SERIES

- 10x eyepiece with pointer, rotatable monocular head
- 4x, 10x & 40xR objectives, movable stage
- Choice: 15w tungsten or cordless LED lighting with charger
- 0.65 N. A. condenser
- Focusing: coarse only or separate coarse & fine
- Hinged at base to adjust angle
- Microscopes by Walter Products, limited 5-year warranty



Catalog #	Description	List Price	Sale
WA2053L	Coarse only focus / Tungsten	\$149 ⁰⁰	\$129⁰⁰
WA2053-RC	Coarse only focus / Cordless LED	\$189 ⁰⁰	\$169⁰⁰
WA2053FCL	Coarse & Fine focus / Tungsten	\$159 ⁰⁰	\$139⁰⁰
WA2053FC-RC	Coarse & Fine focus / Cordless LED	\$199 ⁰⁰	\$179⁰⁰

2057 SERIES

- 360° rotatable monocular head with 10x eyepiece
- Choice: corded LED or cordless LED
- 4x, 10x, & 40x objectives with reverse nosepiece
- 6-hole disc diaphragm, round movable stage
- Coarse only or coaxial coarse & fine focusing
- Microscopes by Walter Products



Catalog #	Description	List Price	Sale
WA2057L	Coarse / Corded LED	\$142 ⁰⁰	\$119⁰⁰
WA2057RC	Coarse/Cordless LED	\$155 ⁰⁰	\$135⁰⁰
WA2057CXL	Coarse & Fine Corded LED	\$159 ⁰⁰	\$139⁰⁰
WA2057CXRC	Coarse & Fine Cordless LED	\$179 ⁰⁰	\$149⁰⁰

3000 SERIES SILVER & GOLD

- Siedentopf binocular head, 30° inclined, 360° rotatable
- Finite Plan DIN objectives: 4x, 10x, 40xR & 100xR *or* Infinity High Contrast Plan objectives: 4x, 10x, 20x, 40xR & 100xR
- Reversed nosepiece
- Coaxial coarse & fine focusing
- Both 10x eyepieces diopter adjustable
- Built-in mechanical stage
- Built-in LED Kohler illumination with field diaphragm
- With light intensity control
- 1.25 Abbe condenser, adjustable aperture, iris diaphragm with graduated scale
- Microscopes by Walter Products



Call to receive our **FREE** Buyer's Guide Catalog

Catalog #	Description	List Price	Sale
WA-3007-L5 <i>Silver Series</i>	Binocular with Finite Plan 4x, 10x, 40xR, & 100xR	\$799 ⁰⁰	\$649⁰⁰
WA-3008-L5 <i>Gold Series</i>	Binocular with Infinity HC 4x, 10x, 20x, 40xR, & 100xR	\$949 ⁰⁰	\$799⁰⁰

3050 SERIES

- Monocular, 10x eyepiece & pointer
- 45° inclined head, rotatable 360°
- 4x, 10x, & 40xR objectives
- Separate or coaxial focusing
- 3 1/2" round movable stage
- Choice: corded LED or cordless LED
- Positive click stops
- Select models include 100xR objective, square mechanical stage, & Abbe
- Microscopes by Walter Products



Catalog #	Description	List Price	Sale
WA3050-LED	Corded LED / 4x, 10x, 40xR Separate Focus / Movable Stage	\$179 ⁰⁰	\$159⁰⁰
WA3050-RC	Cordless LED / 4x, 10x, 40xR Separate Focus / Movable Stage	\$189 ⁰⁰	\$169⁰⁰
WA3050-100-LED	Corded LED / 4x, 10x, 40xR, & 100xR Coaxial Focus / Mechanical Stage	\$389 ⁰⁰	\$289⁰⁰
WA3050-100-RC	Cordless LED / 4x, 10x, 40xR, & 100xR Coaxial Focus / Mechanical Stage	\$399 ⁰⁰	\$299⁰⁰

BMT SERIES

- Monocular, 10x eyepiece with calibrated pointer
- 45° inclined, rotating 360° head
- 4x, 10x, & 40xR objectives
- Reversed nosepiece
- Separate coarse & fine focusing
- Light - choice of corded LED or cordless LED
- Includes variable light intensity control
- 3/4" x 1/2" stage with stage clips
- 0.65 N.A. condenser with 5-hole disc diaphragm
- Rugged cast metal stand construction with stain resistant finish
- Includes cord holder, dust cover & instruction manual
- Limited 5-year warranty
- Microscopes by Walter Products



Catalog #	Description	List Price	Sale
WA-BMT-LED	Corded LED	\$239 ⁰⁰	\$199⁰⁰
WA-BMT-RC	Cordless LED	\$249 ⁰⁰	\$209⁰⁰

Visit our Website



40 SERIES

- 10x eyepiece & pointer, rotatable monocular or binocular head
- 4x, 10x & 40xR objectives or 4x, 10x, 40xR & 100xR
- Choice: corded or cordless LED lighting with charger
- Coaxial coarse & fine focusing, iris diaphragm
- Combination mechanical stage, Abbe, & 100xR on select models
- Microscopes by Walter Products, limited 5-year warranty



Catalog #	Description	List Price	Sale
WA40CXM Monocular	Corded LED / 4x, 10x, 40xR 0.65 N.A. condenser / Plain Stage	\$269 ⁰⁰	\$239⁰⁰
WA40CXMRC Monocular	Cordless LED / 4x, 10x, 40xR 0.65 N.A. condenser / Plain Stage	\$289 ⁰⁰	\$259⁰⁰
WA40CXM100 Monocular	Corded LED / 4x, 10x, 40xR, 100xR 1.25 Abbe condenser / Mechanical Stage	\$369 ⁰⁰	\$319⁰⁰
WA40CXM100RC Monocular	Cordless LED / 4x, 10x, 40xR, 100xR 1.25 Abbe condenser / Mechanical Stage	\$389 ⁰⁰	\$339⁰⁰
WA40CXB 100 Binocular	Corded LED / 4x, 10x, 40xR, 100xR 1.25 Abbe condenser / Mechanical Stage	\$499 ⁰⁰	\$429⁰⁰